Therapy Equipment by Jenx

Ziggy offers a unique experience in positioning.

Ziggy is used for play, support, therapy, rest and imaginative

play. There are many ways to use Ziggy: to promote activity

such as crawling, to promote balance when used for sitting

either astride or to one side, to support a posture by placing

This highly adaptable long, soft and flexible product is square

in cross section and can be bent at a number of points to

create many different shapes. It can be easily transported

too, so if you are away from home and need a corner seat

shape and place on the floor – an instant corner seat.

to support a young child, simply bend Ziggy round into a 'U'

If you want to place your child prone with shoulders raised,

position over Ziggy and bend the product around their hips

but you don't have a wedge, simply place him or her in

Everyone can find a use for Ziggy, the question is:

how many uses for Ziggy can you come up with?

Ziggy behind or in front of a child to help to maintain an

overall body posture and the list goes on.

to keep your child in mid-line.

Jenx Therapy items are bright, easy to use therapy aids that children can quickly start to benefit from without undergoing formal assessment.

The simple design of the Therapy range means they can offer benefits for a number of children within a department rather than after assessment for an individual child. They can help with support or specific therapy treatment such as postural drainage or prone lying, and can be used as part of regular therapy programmes.

The Therapy range includes items that are simple and flexible, whilst providing fun, comfort and achievement for children with all kinds of developmental needs. There are rolls and wedges to climb over, or to rest on or over, and tables and chairs that you can use to take part in all your favourite floor based activities.

Use the nursery table to play games on with a friend, or tilt it into an easel so you can create marvellous works of art and creativity. And Ziggy is a bendy wonder!

Always Encouraging

Therapy equipment can improve your child's outlook on the world and benefit in a variety of ways.











Children gain huge psychological and physical benefits from therapy programmes, using equipment that can aid physical, social and developmental improvements.

Size 1 mm/inch L: 1200mm H: 100mm W: 100mm / L: 471/4 H: 4 W: 4 Size 2 mm/inch L: 1200mm H: 150mm W: 150mm / L: 471/4 H: 53/4 W: 53/4



lexible and *Versatile*

Ziggy can be bent to offer support to bodies limbs and joints in lots of very intricate and specific ways. The shapes are retained under gentle pressure. While its long length means it offers lots of support and can be formed into whole circles of support if necessary.

Ziggy can be used to support almost anyone in any position. It works well with Dreama™ where you wish to support a defined whole body or limb posture.

Wedges, Rolls and Mats

A range of soft, colourful and comfortable supports in various shapes and sizes, which can be used to create play opportunities and for specific therapy treatments.

Can be used for prone and also supine lying and have accessories to enhance postural control and stability.

Easy Care

Plus the easy clean, water repellent, fire

retardant covers make your item safe,

easy to keep clean and pleasant to use.

Wedges: See above diagram for size a,b,c,d area.

Can be used in lying to support both legs at the knees. This tends to increase the stability of the child and helps to prevent legs rolling to the side or the child "windsweeping". T-rolls are particularly effective in reducing windsweep in lying.

Rolls are also a great aid to practice sitting balance and to practice sit to stand, whilst maintaining a neutral or abducted hip position for the child to improve the standing position achieved.

Durable, Colourful Colours

Covers made from durable

abrics (and 5 different colour

options) lengthen the life of your

Always Thoughtful Outstanding features and benefits

Always Configurable Accessories

Size 1 mm/inch a:600mm x b:190mm x c:600mm x d:20mm / a:23% x b:7½ x c:23% x d:3%

Size 2 mm/inch a:500mm x b:150mm x c:500mm x d:20mm / a:19% x b:5% x c:19% x d:%

Size 3 mm/inch a:1200mm x b:350mm x c:600mm x d:20mm / a:471/4 x b:133/4 x c:233/4 x d:34

Size 4 mm/inch a:920mm x b:300mm x c:600mm x d:20mm / a:36¼ x b:11¾ x c:23¾ x d:¾

Size 5 mm/inch a:600mm x b:210mm x c:600mm x d:50mm / a:23% x b:181/4 x c:23% x d:3/4

Mats:

Size

Abduction block to help maintain

a neutral or abducted hip position

in prone lying. Comes in 2 sizes.

Size 1 mm/inch L/D 600mm x 300mm / 23% x 11%

Size 2 mm/inch L/D 900mm x 300mm / 35½ x 11¾

Size 3 mm/inch L/D 1200mm x 300mm / 471/4 x 113/4

Side Walls help prevent

children from rolling off the

Custom made to your specification

These are custom made to the size and shape your child needs.

When you have decided what you need, simply get in touch with us at Jiraffe, and we'll send you out a form which will allow you to specify measurements and other requirements, then we will provide you with a full quotation.

Range of Accessories

Wedges used for prone lying can have

children from rolling off the wedge and

abducted hip position in prone lying.

Wedges used for supine lying can also

have straps, side walls and abduction and

may well be used in conjunction with a roll

to offer a comfortable position for the legs

whilst reducing abdominal cavity pressure,

especially helpful immediately after a meal.

for Wedge

Straps keep a child

centred on the wedge

abduction blocks to maintain a neutral or

straps to keep a child centred on the

wedge, side walls to prevent active

Please bear in mind that custom-made items cannot be returned so do take great care with your measuring.

Always Thoughtful Outstanding features and benefits

Corner Seat by Jenx

The Corner Seat is available in

It is an adjustable floor based seating

offers support in long sitting, which is

children's usual first sitting position.

three sizes to suit children from six

months to approximately 10 years.

system for infants and young children and

The Corner Seat has an adjustable height

back and can be used with or without the

upper back support. An abduction block is

used to help keep your child well positioned,

with their hips neatly back into the seat.



Activities

Can be used with the Jenx Nursery Table to create a great working space for lots of fun activities including working and playing with



child in midline and, if you prefer to have a flat back with wings at the side, there are small triangular cushions available to place at the back to provide this setting.

The Corner Seat is a bright, comfortable and ideal first seat for any child, and can be used with the Jenx Nursery Table. Small touch and close tape straps are used to keep the two products safely connected to each other whilst in use.

The corner shape helps to maintain your

Trunk **Positioning**

Sturdy 'wings' and corner design help to maintain symmetry of your child's spine.

other children across the table

Always Configurable Accessories



For abduction when long sitting, the abduction block is adjustable to ensure the right position for the child.

The Jenx Corner Seat has various accessories to ensure it meets every child's needs.

Simple first seat to

ntroduce children

o floor sitting. Long

sitting, with the legs

stretched out in front

helps to maintain the


length of the muscles

behind the knee.

Corner Seat Abduction block

Corner Seat Triangular Cushions For those who prefer to have a flat back with wings. The triangular cushions can also be used beside the legs to prevent external rotation of the legs.





Nursery Table by Jenx

The Jenx Nursery Table is available in three sizes to suit children from nine months to six years old.

It has a large, plastic play/activity surface and can be used flat, angled or as an easel. The frame is adjustable for height and the Nursery Table can be used to great effect with the Jenx Corner Seat. A short strap attaches the two items in use to prevent a child tipping the table over.

The plastic surface is easy to wipe clean and it is cut out to one side to allow the child to get really close to his or her table-top activities, and has edging around the remainder of surface to reduce the risk of objects sliding off. The Jenx Nursery Table is cheerful and adaptable – a great item for all sorts of activities.

Sturdy and Bright

writing activities.

Useful Cutout Detail

The cut out on the table gives children the opportunity to be surrounded by the work surface, which means that even children who find it difficult to bring their hands well forward will have their arms supported. This in turn helps them to access the activity on offer. Many young children find the cut out

gives them a feeling of security. Children can

work from either side of the table, making

this an ideal system for children to play on

out side and the other on the side with the

together. One child can work on the cut

Sturdy table with large, clear surface area.

The large surface is big enough for a child

to play or work at most activities with ease.

The surface is ideal for using for drawing and

Bright contrasting table edging, which assists

children in identifying their area of work. It pre-

vents items from being knocked off the table

and provides a place for toys to be secured

for dirt to be brushed through for cleaning.

to. Also the space beneath the edging allows

Always Thoughtful Outstanding features and benefits

Height Adjustable The Nursery Table can be

easily adjusted for height. Adjustment is continuous rather than incremental, so every child can have their table at precisely the correct height for them. This feature allows the Nursery Table to be used by children using a Corner Seat or other floor seat, or by children who are able to sit unaided. This gives rise to endless play

opportunities with

friends.

Corner Seat Tape Straps for attaching the Jenx Nursery Table, so your child can enjoy all sorts of fun activities.

> Visit **jenx.com** for additional product details.

Advanced Side Lying Board by Jenx

Side Lying Boards offer supported comfort to children while resting or during play activities. Side lying is a great position in which to recharge and also for joining in with classroom activities.

The Jenx Advanced Side Lying Board is not only comfortable, but both the lying and back support surface can be adjusted to prevent children rolling forward. An optional support pack is also available which comprises head support and two leg supports to offer support for the legs to help to maintain a neutral hip position.

Larger sizes of the Jenx Advanced Side Lying Board come with a mobile plinth to raise the lying surface to waist height, which means ease of hoisting, less manual handling for carers and also means children can interact with peers more easily.

Jenx Advanced Side Lying Boards offer great support for children who have poor postural control, low muscle tone or have a tendency to curve to the side. In side lying the back is fully supported preventing the spine from curving. The generous upholstery means that children can be safely and comfortably supported during rest or activities. Sturdy, adjustable straps come as standard.

Always Thoughtful Outstanding features and benefits

Comfortable, Adjustable Back Support

Full length support for a child means they feel safe and can relax. The back can be angled to help to prevent your child rolling forwards.

Comfortable Adjustable

Lying Surface The well-upholstered surface is comfortable to lie on and wide enough for your child to feel secure. The angle of the surface can be altered to slope your child slightly backwards to prevent him or her rolling forward.



For added security, wide straps are used to keep your child well positioned and to add to his or her safety. The straps will normally be used around the chest

and hips, but each child


may need support in a slightly different way.

Jenx Advanced Side Lying Board Technical Specification

Side Lying Boards	Units	SL2P Size 2	SL3P Size 3	SL4P Size 4
Age Range (Approx)	Months/Years	9mths - 4yrs	3yrs - 9yrs	8yrs - 18yrs
Base Size	mm/inch	950 x 410 / 37½ x 16	1250 x 480 / 49¼ x 19	1860 x 680 / 731/4 x 263/4
Max User Weight	Kg/lbs	80 / 1761/4	80 / 1761/4	80 / 1761/4
Bench Length	mm/inch	1335 / 52½	1640 / 64½	1935 / 76
Bench Width	mm/inch	555 / 22	650 / 281/2	740 / 29
Bench Angle (Supine/Prone)	Degrees	15°/0°	15°/0°	15°/0°
Back Angle (Supine/Prone)	Degrees	25°/90°	25°/90°	25°/90°
Back Height from Bench	mm/inch	330 / 13	350 / 13¾	400 / 15¾
Floor to Bench Height	mm/inch	590 / 231/4	625 / 24½	760 / 30

Bright and Colourful


The Therapy range from Jenx has the ability to brighten everyone's day with a range of 5 sizzling shades to choose from:

















Always here

To arrange a demonstration or assessment with the Jenx Therapy range, just get in touch!

My name is:

Please ring me on:

Or email me at:







Therapy Equipment by Jenx A colourful, simple and flexible range of products to

help with support or specific therapy treatment



Always Supportive...

Jenx Limited, Wardsend Road, Sheffield S6 1RQ, United Kingdom tel: +44 (0) 114 285 3376 email: export@jenx.com web:www.jenx.com



Jenx Limited, Wardsend Road, Sheffield S6 1RQ, United Kingdom tel: +44 (0) 114 285 3376 email: export@jenx.com www.jenx.com



